



# DreamMapper 2.4

## Mobile App User Guide

June 8, 2016

## Table of Contents

<b>Sleep Screen</b> .....	<b>1</b>
How to navigate .....	2
Changing Data Categories .....	5
See Additional Details for the Day and Category .....	6
<b>Other Menu Choices</b> .....	<b>7</b>
<b>Feed Screen</b> .....	<b>8</b>
<b>Learn Screen</b> .....	<b>9</b>
Apnea .....	9
Equipment.....	10
Troubleshooting.....	11
<b>Goals Screen</b> .....	<b>12</b>
<b>Coaching Screen</b> .....	<b>13</b>
<b>Reminders Screen</b> .....	<b>14</b>
<b>Settings Screen</b> .....	<b>15</b>
Account .....	15
Equipment.....	15
Contact Preferences.....	15
<b>Help Screen</b> .....	<b>16</b>
FAQ.....	16
Support .....	17
About.....	18



Welcome!

Welcome to DreamMapper<sup>1</sup>. Here are a few tips to familiarize you with this exciting app.

## Sleep Screen

When you first log on to DreamMapper the “Sleep” screen is displayed and it defaults to the current date. Data is presented in one of three categories: AHI, Usage, and Mask Fit. Usage is the default category. The graph below shows the values for the selected category for a 7-day period including the displayed day. The graph highlights the currently-selected day on the graph in dark blue.



<sup>1</sup> DreamMapper is the next generation of the SleepMapper app from Philips



## How to navigate

### Changing the Date: Tap the chevrons

To go backward or forward in time, tap the chevrons (“<”, “>”) located beside the date. As you change the date, the highlighted bar in the graph and the details at the bottom of the screen reflects the selected date.



### Changing the Date: Interacting with the graph

The graph responds to your touch, so that you can easily see your data for a different day or even a different week.



When a specific day of the graph is tapped, the focus of the currently selected Category will change to the tapped day.





Swipe left-to-right on the graph to see the preceding week, or right-to-left to see the following week. The displayed date range confirms the new time period that's being presented by the graph.






## Changing Data Categories

To change the category and graph: tap the desired category


The selected Category will always be in **bold**, and the default is Usage. To view your AHI or Mask Fit information for a different category tap on the label or value for the category. The pointer at the top of the graph section confirms your selection.





[See Additional Details for the Day and Category](#)

The bottom portion of the Sleep screen will show additional details pertaining to the Category you currently have selected in context of that day.








## Other Menu Choices

To navigate to the other menu items tap the menu icon


To see the **Menu** in the mobile apps tap the icon on the top-left corner of your screen. When tapped the menu will appear, revealing other areas of the application: Feed, Learn, Goals, Coaching, Reminders, Settings, and Help.





## Feed Screen

The Feed section contains reminders, alerts, and recommendations sent to your account.






## Learn Screen

The Learn section contains links to other helpful information related to sleep apnea, equipment usage, and troubleshooting.


### Apnea

This page contains links to several videos addressing the causes and treatment of sleep apnea.



## Learn


[Apnea](#)[Equipment](#)[Troubleshoot](#)



### What is OSA

Duration: 1:48


This video provides an overview on obstructive sleep apnea.



### Sherry and Ashley from the Biggest Loser discuss their OSA stories

Duration: 6:10


Ashley and Sherry from the Biggest Loser discuss their experience with obstructive sleep apnea.



### Overcoming Issues with Therapy

Duration: 1:54

How to overcome possible issues that may occur from PAP therapy treatment for obstructive sleep apnea.



### Living with therapy

Duration: 3:33




## Equipment

This page contains links to videos and online guides related to the use of therapy equipment and DreamMapper.

 **Learn**

[Apnea](#) [Equipment](#) [Troubleshoot](#)

**Guides**



**DreamStation simple device operation**  
Duration: 1:48  
How to operate DreamStation simple device

[Mask Guide](#) >

[Replacement Guide](#) >

[First Night Guide](#) >

[Welcome to DreamMapper User Guide](#) >

**Videos**



**Wisp pediatric fitting**  
Duration: 2:35  
Video providing fitting instructions for the Wisp pediatric mask



## Troubleshooting

This page contains links to videos on the subjects of mask issues, device pressure, equipment issues, and side effects.

Learn	
Apnea	Troubleshoot
Mask Issues	
Correcting Mask Leaks	>
Getting Rid Of Pressure Points	>
General Discomfort	>
Device Pressure Issues	
Mask and Tube Settings	>
Adjusting Ramp Settings	>
Adjusting Ramp Start Pressure	>
Equipment Issues	
Tubing Management	>
Adjusting Humidifier	>
Side Effects	
Nasal Issues	>



## Goals Screen

The Goals screen displays several categories of monthly and personal goals.



Click on the ">" symbol to the right of each category to view its definition and to set your personal goals.



## Coaching Screen

This section contains links to mini-questionnaires and videos pertaining to motivation, heart risk, and your feelings about sleep apnea and sleep therapy.


Coaching	
Choose a topic	
Motivation	>
Heart Risk	>
Feelings	>

Click on the ">" symbol to the right of each topic to view the mini-questionnaire and video.



## Reminders Screen

The Reminders screen is where you set the time intervals at which you will be reminded to clean your mask, humidifier, tubing, and device filters. You can also schedule follow-up reminders. The reminders you schedule here will appear on the Feed screen.



Reminders

Cleaning Reminder

**Mask reminder**

☒

We recommend that you clean your mask every day.

Every Other Day

**Humidifier reminder**

☒

We recommend that you clean your humidifier every week.

Weekly

**Tubing reminder**

☒

We recommend that you clean your tubing every other day.

Daily

**Device Filters reminder**

☒

We recommend you clean your filters every month.

Daily

Follow-Up Reminder

**Follow-Up Reminder**


☐

Most insurance providers require you to follow-up after using therapy for a certain number of days. Check with your homecare provider for more details.





To “turn on” a reminder, click on its  button. The button changes to this: 

Next, select a time interval from the drop-down list: 

For the Follow-up Reminder, select or enter a date in the space provided.

## Settings Screen

The Settings section contains information related to your account, preferences, mask, and data connections.

### Account

From here you can change your Username or password. You can also log out of DreamMapper.



### Equipment

The options in this section allow you to identify your mask, therapy device, and if you want to use Bluetooth communications.

### Contact Preferences

Use the options in this section to specify how DreamMapper should contact you with important information about your therapy data.




## Help Screen

The Help section contains general information about DreamMapper, terminology regarding sleep apnea and sleep therapy, and Philips contact information.

### FAQ


Some of this information may be new or unfamiliar to you. In the FAQ section you will find definitions for all the fields displayed within DreamMapper.






## Support

If you still have difficulty using DreamMapper, please look under Help-Support to find Philips' contact information. We want you to love DreamMapper as much as we do.




# Help

[FAQ](#)[Support](#)[About](#)



## Philips Respironics


Respironics Inc.  
1001 Murry Ridge Lane  
Murrysville, PA 15668 USA



Respironics Deutschland  
Gewerbestrasse 17  
8211 Herrsching, Germany

EC

REP



1 (855) 699 - 6276  
1 (724) 387 - 5036  
[dreammapper.support@philips.com](mailto:dreammapper.support@philips.com)



## About

This section provides the intended use of DreamMapper, applicable cautions and warnings, symbol definitions and additional Philips contact information.

The screenshot shows a mobile application interface for 'DreamMapper'. At the top is a blue header with a white menu icon and the word 'Help'. Below the header is a navigation bar with three tabs: 'FAQ', 'Support', and 'About', with 'About' being the active tab. The main content area is titled 'About DreamMapper' and contains the following text:

The DreamMapper software allows patients on supported Philips Respironics PAP therapy devices to view their compliance data through their smartphone (Mobile Application) or Internet access (Web Application). DreamMapper provides compliance, coaching, and motivational feedback to encourage habitual use of Philips Respironics PAP therapy devices. It also allows data to be transmitted between supported Philips Respironics PAP therapy devices and the healthcare provider. DreamMapper is highly secure to help meet requirements for HIPAA compliance and privacy laws. This includes access security as well as the security and privacy of transmitted data.

**Intended Use**

DreamMapper is intended to support patients by displaying therapy data recorded for patients who are prescribed compatible therapy devices in accordance with the intended use of those therapy devices. DreamMapper provides remote viewing access to patient data via their smartphone (Mobile Application) or Internet access (Web Application) to understand therapy usage, provide encouragement, and coaching towards sufficient usage of supported Philips PAP therapy devices. The DreamMapper software does not perform automatic scoring or diagnosis. The viewable data is to help the patient understand their therapy usage and not to diagnose any condition of the patient.

**Password Policy**

Passwords are required to be between 8 and 20 characters in length.

Passwords must have 3 of the following:

- Uppercase characters
- Lowercase characters
- Numerical characters
- Special characters (e.g., #5&@)

Passwords cannot contain your email address or the word "password".



## Help – About continued:

### Cautions

U.S. Federal law restricts this device to sale by, or on the order of, a licensed physician.

### Warnings

DreamMapper should not be used to diagnose any condition of the patient. Associated therapy device manuals should be read and understood prior to viewing compliance data. Use DreamMapper with Philips Respironics supported therapy devices only. DreamMapper is not compatible with therapy devices not supported by Philips Respironics. If you are not sure if your therapy device is supported by DreamMapper, contact Philips Respironics Product Support for information on the supported devices.



Read and understand the associated therapy device manuals prior to viewing compliance data.

Respironics Inc.  
1001 Murry Ridge Lane  
Murrysville, PA 15668 USA



Respironics Deutschland  
Gewerbestrasse 17  
82111 Herrsching, Germany



Toll-free Telephone : 1 (855) 699 - 6276  
Telephone : 1 (724) 387 - 5036  
Email : dreammapper.support@philips.com

### App Version

2.4 (Alpha May 20) (64)

1119424 R07  
RPM 5/11/2016



1101478, 1101479, 1101679, 1128262, 1128264,  
1128265

[License Agreement](#)



[Privacy Notice](#)





Respironics, Inc.  
1001 Murry Ridge Lane  
Murrysville, PA 15668 USA



Respironics Deutschland  
Gewerbestrasse 17  
82211 Herrsching, Germany

EC REP



REF

1101478, 1101479, 1128264,  
1128265  
1125925 R01  
RPM 06/08/2016