

# CoughAssist T70

## Prescription

Patient name \_\_\_\_\_

Date of birth \_\_\_\_\_

Date \_\_\_\_\_

### Settings

Mode:  manual  auto

Interface:  mask

Cough-Trak:  on  off

mouthpiece

trach adaptor

Inspiratory pressure (range) \_\_\_\_\_ cm H<sub>2</sub>O

Inspiratory time \_\_\_\_\_ secs.

Expiratory pressure (range) \_\_\_\_\_ cm H<sub>2</sub>O

Expiratory time \_\_\_\_\_ secs.

Frequency of treatment \_\_\_\_\_

Pulse oximetry SpO<sub>2</sub> ≤ \_\_\_\_\_ %  NA

Titrate to achieve inspiratory pressures of 35-40 cm H<sub>2</sub>O and expiratory pressures of -35-45 cm H<sub>2</sub>O or to achieve an effective cough.

### Additional orders

### Duration of use

lifetime

other \_\_\_\_\_

### Physician information

Name (please print) \_\_\_\_\_

Signature \_\_\_\_\_

Telephone \_\_\_\_\_