



Therapeutic support

Facilitating developmentally supportive care practice

PHILIPS

Children's  Medical Ventures

sense and simplicity



Safety is a critical element in any healthcare situation but it is especially important when providing care to fragile and vulnerable premature infants. To assist clinicians in locating the tiny veins in the premature infant's hand, Philips Children's Medical Ventures offers the Wee Sight Transilluminator. This device helps to prevent multiple insertions.

In the developing fetus, systems develop in a set order:

- Tactile
- Vestibular
- Gustatory
- Olfactory
- Auditory
- Visual

When the premature infant enters the world, many of these systems are not yet developed. It is important to provide protection to enable the continued development of these systems, outside of the uterus. For example, the visual system is the very last system to develop at around 32 to 34 weeks gestation. For any infant born prior to 35 weeks, it is critical to protect their eyes so that the visual system can continue to develop. To help in protecting the infant's eyes, consider lowering the ambient light level and providing eye protection, especially during the following procedures:

- Pre and post ROP examinations
- Line insertions
- Any other procedure or exam requiring bright procedure lights

Philips Children's Medical Ventures offers a variety of eye protection to meet your needs.



Wee Sight Transilluminator

Answers in the palm of your hand

Wee Sight Transilluminator can help to accurately locate the tiny veins in a neonate's small limbs for improved insertion of an IV. The strong LED light performs as well as larger transilluminators, but does not emit heat – making it safer even for the most delicate skin.

The Wee Sight:

- Uses two AA alkaline batteries (included)
- Is a Class 2 LED product
- Has an output light wavelength of 629 nm

Product description	Quantity	Part number
Wee Sight Transilluminator	1	1017920
Heart-patterned lanyard	1	1026486
Solid purple lanyard	1	1017291



**Philips Healthcare is part of
Royal Philips Electronics**

How to reach us

www.philips.com/healthcare

Fax: +31 40 27 64 887

Asia

+852 2821 5888

Europe, Middle East, Africa

+31 40 27 87246

Latin America

+55 11 2125 0764

North America

+1 425 487 7000

+1 800 285 5585 (toll free, US only)

Philips Children's Medical Ventures

191 Wyngate Drive

Monroeville, PA 15146

+1 412 380 8800

+1 800 345 6443 (toll free, US only)

Customer Service

+1 724 387 4000

+1 800 345 6443 (toll free, US only)

Respironics Europe, Middle East, Africa

+33 1 47 52 30 00

Respironics Asia Pacific

+852 3194 2280

www.philips.com/childmed



© 2009 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication. Respironics, Bendy, SnuggleUp, Frederick T. Frog, Cozycare, Prone Plus, Gel-E Donut, Wedgie, Squishon, IncuCover, Tucker Sling, Wee Thumbie, Wee Soothie, NICU Soothie, Super Soothie, Soothie, Sweet-Ease, Snoedel, WeePee, Pee-On, Sali-Wipes, Tiny Tub, Heel Snuggler, Wee Sight, WeeSpecs, Global Car2e, Wee Focus and Premie for a Day are trademarks of Respironics, Inc. and its affiliates. Air-Shields is a trademark of Dräger and Care Plus, Giraffe, and OmniBed are trademarks of GE Medical.

RIT CHMV 009 E * JAN 2009 bdc10994