

Philips DreamWear  
Nasal mask

Fit pack — frame and 4 cushions  
INTL

HH1111



## DreamWear®

Nasal mask with mask cushion underneath the nose

The innovative design\* of DreamWear offers you the freedom and comfort to sleep in the way you would like.

### **Freedom of movement**

- Innovative mask design for more freedom in sleep therapy
- More comfort with DreamWear
- For a more natural sleep
- Soft, flexible mask frame

### **No red marks on the bridge of the nose**

- Cushion under the nose
- Slim headgear design

### **Unobstructed view**

- Hose connection on the head
- Read without restrictions or watch TV

**PHILIPS**

Nasal mask  
Fit pack — frame and 4 cushions INTL

# Highlights

## Cushion under the nose



The mask cushion lies gently under the nose, which causes no pressure points on the bridge of the nose and prevents unpleasant skin irritations on the nostrils.

## A more natural sleep



DreamWear combines the advantages of minimal-contact masks and nasal masks, provides high comfort and can contribute to successful sleep therapy.

## Hose connection on the head




The exhaled air is not discharged directly at the nose, but at the top of the head. This hose positioning allows more freedom of movement during sleep. Users say they no longer need to sleep in a specific position with DreamWear, but have more freedom.\*

## Read or watch TV



Thanks to the open design, DreamWear offers a very large field of view, so you can easily wear glasses and read, watch TV or use a tablet PC in bed before you go to sleep.

## Soft, flexible mask frame



The flexible material through which the air flows gently rests on the skin, preventing pressure points and providing a comfortable fit.

## Innovative mask design

The innovative design\* of DreamWear prevents pressure sores and skin irritation on the nostrils and on the nose bridge, and allows the user more freedom of movement through the hose connection on the head, as well as a clear field of view.

## Slim headgear design



DreamWear unites comfort and aesthetics — the headgear is unobtrusive and sleek in colour, making DreamWear as discreet and unobtrusive as possible.

## More comfort with DreamWear

Users said DreamWear was more comfortable to wear than their previously prescribed mask, as it is barely noticeable.

HH1111/01

# Specifications

## Freedom of movement

- Clear field of vision: Individual sleep position

## Innovative design

- Comfort and performance: Reduced skin irritation

## Service

- Warranty: 90 days

## Intent of use

- To be used on patients >30 kg (>66 lbs): with CPAP prescription

## Specifications

- Operating pressure: 4–20 cmH<sub>2</sub>O
- Sound pressure level: 18 dBA
- No. of parts: Five
- Headgear material: Urethane foam/nylon/Spandex
- Frame material: Polycarbonate
- Cleaning instructions: mask: Hand wash with mild detergent, Rinse thoroughly, air dry
- Cleaning: headgear, tubing: Hand wash with mild detergent, Rinse thoroughly, hang to dry
- Frame/Elbow/Swivel material: Polycarbonate
- Fabric sleeve material: Polyester/Spandex



Issue date 2018-05-10

Version: 6.3.1

12 NC: 0000 011 16700  
UPC: 6 06959 02311 6

© 2018 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* Source: Philips User Survey 2015

\* The costs of DreamWear are private and are usually not reimbursed by statutory health insurance.

\* The scope of supply includes a standard mask frame (size M) and four cushions for all sizes (S, M, MW, L). The sleep therapy device is not included.