



Philips DreamWear
Full-face magnet clips

Full-face replacement part
with 1 set of magnetic clips



HH1132

Like wearing nothing at all.*

With the freedom and comfort of DreamWear Full Face, you can sleep how you want.

Magnetic clips

- The headgear clips and mask cushion contain magnets

Easy to use

- Magnet clips for headgear attachment

PHILIPS

Full-face magnet clips
Full-face replacement part with 1 set of magnetic clips

Specifications

Warning: Magnetic interference

- Beware with pacemakers: Read IFUs before use.

Product details

- Warranty: 90 days

HH1132/00

Highlights

Easy to use



The magnet clips make attaching the headgear to the frame easy and fast!

Magnetic clips

Contact your healthcare professional before you use this mask. Some medical devices may be affected by magnetic fields. The magnetic clips in this mask should be kept at least 50 mm (2 in.) away from any active medical device with special attention to implanted devices, such as pacemakers, defibrillators and cochlear implants. Read in-pack IFUs for more information.



Issue date 2018-06-22

Version: 2.2.1

12 NC: 0000 011 33454
UPC: 6 06959 41466 2

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

** Data analysis of Aug 2017 Patient Preference trial where n = 85 and prescribed masks include Resmed Airfit F10 and F20, Resprionics Amara View and F&P Simplus; based on 10 days of use.