



Philips DreamWear
Minimal contact nasal mask

Pack
and 4 cushion sizes
3 frame sizes

HH1121

Like wearing nothing at all*

With the freedom and comfort of DreamWear, you can sleep how you want *2015

Philips User Preference Questionnaire

Innovative design

- Air tube connection at the top of the head
- Under-the-nose cushion
- Supports your sleep therapy
- Four-part modular design

Open field of vision

- Open design offers an extremely wide field of vision

Comfortable design

- Soft, flexible frame
- Slim-line headgear
- Soft fabric wraps



PHILIPS

Minimal contact nasal mask
Pack and 4 cushion sizes, 3 frame sizes

HH1121/01

Highlights

Freedom of movement



By moving the tube location to the top of the head, DreamWear offers more freedom of movement during the night. Users stated they didn't need to choose a sleep position when wearing DreamWear.*

Four-part modular design



Modular design for quick assembly and cleaning, with interchangeable cushions available in four sizes.

Slim-line headgear



DreamWear combines comfort with aesthetics: Not only will you barely feel the mask, it is also subtle and unobtrusive.

Soft fabric wraps



Provides a soft interface to the cheeks

Soft, flexible frame



The flexible material through which the air flows is applied gently to the skin and is pleasant to wear.

Supports your sleep therapy



DreamWear combines the advantages of minimal contact masks and nasal masks. It is comfortable and can help a successful start to sleep therapy.

Under-the-nose cushion



Prevents red marks, discomfort or irritation in the nostrils or on the nose bridge. Users rated DreamWear more comfortable than their prescribed mask.

Unobstructed view



It's easy to wear glasses, read, watch TV, and use a computer or tablet before falling asleep.

Specifications

Service

- Warranty: 90-days

Innovative design

- Open field of vision: Wear glasses, read or watch TV

Comfort

- Under-the-nose cushion: Prevents red marks in nostrils

Intent of use

- To be used on patients >30kg: with CPAP prescription

Specifications

- Operating pressure: 4 -20 cmH₂O
- Sound pressure level: 18 dBA
- N° of parts: Five
- Cushion material: Silicone
- Headgear material: Urethane foam/nylon/Spandex
- Frame material: Polycarbonate
- Cleaning instructions: mask: Hand wash with mild detergent, Rinse thoroughly, air dry
- Cleaning: headgear, tubing: Hand wash with mild detergent, Rinse thoroughly, hang to dry
- Frame/Elbow/Swivel material: Polycarbonate
- Fabric sleeves material: Polyester/Spandex



Issue date 2018-04-03

Version: 3.1.1

12 NC: 0000 011 16695
UPC: 6 06959 02298 0

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* Source: Philips User Survey 2015

* The purchase of this medical device does not give the right to
obligatory Social Security reimbursement