

Philips Wisp  
Magnetic clips

2 per pack

HH1027



## Clips make the headgear easy to adjust

Click-fit headgear clips with adjustable tabs. Allows quick and easy fitting and adjustment.

### Easy to use

- The headgear clips and mask cushion contain magnets

### Design

- Click-fit headgear clips with adjustable tabs

**PHILIPS**

Magnetic clips  
2 per pack

# Specifications

## Warning: Magnetic interference

- Beware with pacemakers: Read IFUs before use.

## Service

- Warranty: 90 days

## Design

- Easy fitting: Self-adjusting

## Specifications

- Headgear clips material: Acetal copolymer and neodymium magnet

HH1027/00

# Highlights

## Easy to adjust



A minimal part count and click-fit headgear clips make it easy to use and maintain.

## Magnetic clips

Contact your healthcare professional before you use this mask. Some medical devices may be affected by magnetic fields. The magnetic clips in this mask should be kept at least 50 mm (2 in.) away from any active medical device with special attention to implanted devices, such as pacemakers, defibrillators and cochlear implants. Read in-pack IFUs for more information.



Issue date 2018-04-24

Version: 4.1.1

12 NC: 0000 011 20434  
UPC: 6 06959 03969 8

© 2018 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)