Philips LiquiCell nasal CPAP cushion

Reduces skin irritation Soft and comfortable 30-pack

HH1401

# For a comfortable mask fit

### Ultra-thin and liquid-filled cushions

A soft and comfortable liquid cell that helps reduce skin irritation. When applied across the nasal bridge area, LiquiCell helps to reduce skin irritation that may sometimes occur when using positive airway pressure (PAP) therapy.

Compact design • Contains 30 CPAP cushions Comfortable • LiquiCell acts a pad between the mask cushion and the skin. Medical technology • Protects your skin



## Specifications

Design Soft and comfortable material: Ultra-thin and liquid cushions

### HH1401/00

## Highlights

 $\begin{array}{l} \label{eq:comfortable liquid cell} \mbox{LiquiCell}^{\ensuremath{\mathbb{R}}} \ \mbox{is a medical technology that helps prevent} \end{array}$ skin tissues from being stretched or torn. The ultrathin cells contain a liquid that is in constant motion allowing your body to glide or float with an object—both reducing shear and friction of the skin and improving blood perfusion to the tissues.

Helps reduce skin irritation The LiquiCell Nasal CPAP Cushions are an option for those who are prone to red marks or bumps forming on the bridge of the nose from mask irritation.

### Multipack

LiquiCell is sold as a pack of 30 single-use disposable cushions.

PHI	LIPS
Æ	
E	
	+

Issue date 2017-11-29 © 2017 Koninklijke Philips N.V. All Rights reserved. Version: 4.1.1 Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners. 12 NC: 0000 011 09097 UPC: 6 06959 01271 4 www.philips.com