



Philips DreamWear
Full-face headgear

Full-face replacement part
Standard headgear



HH1131

Like wearing nothing at all.*

With the freedom and comfort of DreamWear Full Face, you can sleep how you want.

Magnetic clips

- The headgear clips and mask cushion contain magnets

Comfort and ease of use

- Soft fabric headgear with magnet clips

PHILIPS

Full-face headgear
Full-face replacement part Standard headgear

Specifications

Warning: Magnetic interference

- Beware with pacemakers: Read IFUs before use.

Product details

- Materials: Nylon, Polyurethane foam, Spandex
- Warranty: 90 days

Replacement

- Inspect daily for wear: Replace at first sign of wear

HH1131/00

Highlights

Soft and comfortable



The soft fabric headgear is comfortable to wear at night while the magnet clips make it easy to adjust for the right fit.

Magnetic clips

Contact your healthcare professional before you use this mask. Some medical devices may be affected by magnetic fields. The magnetic clips in this mask should be kept at least 50 mm (2 in.) away from any active medical device with special attention to implanted devices, such as pacemakers, defibrillators and cochlear implants. Read in-pack IFUs for more information.



Issue date 2018-06-22

Version: 2.2.1

12 NC: 0000 011 33450
UPC: 6 06959 41465 5

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

** Data analysis of Aug 2017 Patient Preference trial where n = 85 and prescribed masks include Resmed Airfit F10 and F20, Respronics Amara View and F&P Simplus; based on 10 days of use.