



Philips SideStream
Adult face mask

Designed for adults

Replace after 6 months
Use with SideStream Nebulizers
Latex and DEHP free

HH1474/00

Easy and reliable

The SideStream mask is used with SideStream high-efficiency nebulizers. This reusable mask is contoured to comfortably fit the face, and the adjustable elastic strap provides a snug fit. This adult mask must only be hand-washed.

Comfortable fit

- Contours to the face

Designed for comfort

- Helps reduce aerosol leaks around the eyes*

Proven SideStream technology

- For specific use with SideStream nebulizers

PHILIPS

Adult face mask

Designed for adults Replace after 6 months, Use with SideStream Nebulizers, Latex and DEHP free

Specifications

Maintenance

- Replace: After 6 months

Product details

- Materials: Latex and DEHP free
- For use with: SideStream Nebulizers
- Intended Use: Designed for adults

HH1474/00

Highlights

Comfortable fit



The soft and flexible material, combined with an elastic strap helps to ensure a comfortable and snug fit.

Designed for comfort



The SideStream mask incorporates a vent design that directs medication away from the eyes.

Fast efficient drug delivery



Our SideStream uses diamond jet technology, providing consistent drug delivery, time after time.



Issue date 2018-07-27

Version: 1.0.1

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

** Data presented by Respironics at the International Society for Aerosols in Medicine 16th International Congress, June 2007, Tours, France.