

Philips Threshold
PEP

Positive expiratory pressure
For airway clearance
Helps remove lung secretions

HH1333



Pressure that's good for you

With a flow-dependent, one-way valve, threshold devices maintain a constant pressure, regardless of the patients airflow.

Convenient design

- Reliable and easy to clean
- Works in any position

Mouthpiece/mask option

- Mouthpiece/mask option

Intended use

- Intended use of PEP

PHILIPS

PEP

Positive expiratory pressure For airway clearance, Helps remove lung secretions

Specifications

Product details

- Warranty: Replace after 1 year
- Intended Use: As directed by your physician
- Cleaning: Wash in warm soapy water after each

- use. Flush all parts with clean water. Shake out excess water and air dry. DO NOT BOIL OR HEAT.
- Includes: training diary

HH1333/00

Highlights

Mouthpiece/mask option

Mouthpiece/mask option for flexibility

Reliable and easy to clean

Threshold is made of durable, high-impact acrylic.

Works in any position

Effective therapy in any position

Intended use

Threshold PEP expiratory pressure trainer for can be used for airway clearance, bronchiole hygiene and chest physical therapy.



Issue date 2018-09-19

Version: 1.0.1

12 NC: 8670 001 46075
EAN: 03 83730 73501 6

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com