



Philips DreamWear
Full face mask

Full Face Pack

with Small cushion
and Small and Medium Frame



HH1126

Like wearing nothing at all.*

With the freedom and comfort of DreamWear Full Face, you can sleep how you want.

Magnetic clips

- The headgear clips and mask cushion contain magnets

Innovative design

- Freedom of movement
- A modular design for more options

Open field of vision

- Open design offers a wide field of vision

Comfort and ease of use

- Easy to use
- Soft, flexible frame
- Under the nose full face cushion



PHILIPS

Highlights

Freedom of movement



By moving the tube location to the top of the head, DreamWear Full Face offers more freedom of movement during the night. Users stated they didn't need to choose a sleeping position when wearing DreamWear Full Face.*

No more nose bridge red marks



The under the nose innovative design of the cushion means no more red marks, discomfort, and irritation on the nose bridge.**

Unobstructed view



Open design makes it easy to wear glasses, read, watch TV, and use a computer or tablet before falling asleep.**

Like wearing nothing at all



The soft, flexible frame is comfortable to wear at night. Users have stated it feels like wearing nothing at all.*

Freedom of choice



The modular design allows users to choose the cushion that works best for them, all on one frame.

Easy to use



Only a few parts for easy assembly, disassembly and cleaning.

Magnetic clips

Contact your healthcare professional before you use this mask. Some medical devices may be affected by magnetic fields. The magnetic clips in this mask should be kept at least 2 in. (50mm) away from any active medical device with special attention to implanted devices such as pacemakers, defibrillators and cochlear implants. Read in-pack IFU's for more information.

Specifications

Warning: Magnetic interference

- Beware with pacemakers: Read IFU's before use.

Product details

- Materials: Nylon, Polyurethane foam, Polycarbonate, Polyester, Silicone, Spandex
- Warranty: 90 days
- Operating pressure: 4 to 30 cmH2O
- Sound level: 17 dBA
- Number of parts: 5

Cleaning: Mask

- Hand wash with mild detergent: Rinse thoroughly, air dry¹

Cleaning: Headgear, tubing

- Hand wash with mild detergent: Rinse thoroughly, hang to dry¹

Replacement

- Inspect daily for wear: Replace at first sign of wear

Intent of use

- To be used on patients >30kg: with CPAP prescription



Issue date 2018-05-30

Version: 3.2.1

12 NC: 0000 011 33375
UPC: 6 06959 04627 6

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

** Data analysis of Aug 2017 Patient Preference trial where n = 85 and prescribed masks include Resmed Airfit F10 and F20, Respiration Amara View and F&P Simplus; based on 10 days of use.

** Design claim - the design of the mask is such that it does not contact the nasal bridge.

¹ Please refer to Instructions For Use for more detailed cleaning instructions.