

Philips DreamStation Go  
Overnight battery

Maintain therapy off the grid

Compact  
Lightweight, 696 g (1.53 lbs)

HH1450



## CPAP therapy even "off the grid"

For many sleep apnoea patients, travelling off the power grid can create problems with maintaining therapy. Philips developed the DreamStation Go overnight battery.

### **Compact and lightweight**

- Compact and lightweight

# PHILIPS

## Overnight battery

Maintain therapy off the grid Compact, Lightweight, 696 g (1.53 lbs)

# Specifications

### Compatibility

- PAP System: DreamStation Go CPAP, AutoPAP

### Service

- Warranty: 1 year

### Technical information

- Weight: 696 g (1.53 lbs)
- Output voltage: 18-24.6 VDC
- Battery technology: Lithium Ions

- Input Voltage Range: 100-240 VAC, 50/60 Hz, 2.0-1.0 A
- Output power: 50 W
- Run time: > 8 hours
- Minimum life cycle: 500 cycles
- Recharge time: < 5 hours

### Part Number

- Battery Pack: 1133281

HH1450/00

# Highlights

### Compact and lightweight



If you plan on travelling off the power grid, DreamStation Go's overnight battery can help you stay committed to your therapy. The DreamStation Go set at 10 cm of pressure displayed an average run time of 13 hours\* when running on the DreamStation Go overnight battery. Note: Battery run times may vary due to higher CPAP pressure settings, excessive mask leakage, altitude, temperature, charging devices with USB ports and other factors. \*PAP run-time test conditions: CPAP mode, pressure 10 cm H<sub>2</sub>O, 12 mm tubing, leak 37 LPM, 23°C room at an elevation of 503 m (1650 ft).



Issue date 2018-08-31

Version: 5.0.2

12 NC: 0000 011 33281  
UPC: 6 06959 41479 2

© 2018 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)