

# PHILIPS

## PEP

### Threshold

Positive Expiratory Pressure device

For airway clearance

Helps remove lung secretions

HH1333



## Pressure that's good for you

With a flow-dependent, one-way valve, threshold devices maintain a constant pressure, regardless of the patient's airflow.

#### Convenient design

- Reliable and easy to clean
- Works in any position

#### Mouthpiece/mask option

- Mouthpiece/mask option

#### Intended use

- Intended use of PEP

# Highlights

**Mouthpiece/mask option**

Mouthpiece/mask option for flexibility

**Reliable and easy to clean**

Threshold is made of durable, high-impact acrylic.

**Works in any position**

Effective therapy in any position

**Intended use**

Threshold PEP expiratory pressure trainer for can be used for airway clearance, bronchiole hygiene and chest physical therapy.

# Specifications

**Product details**

**Warranty:** Replace after 1 year

**Intended Use:** As directed by your physician

**Cleaning:** Wash in warm soapy water after each use. Flush all parts with clean water.

Shake out excess water and air dry. **DO NOT BOIL OR HEAT.**

**Includes:** training diary

