



DreamStation ordering guide



Ordering convenience

You now have the ability to order your award winning DreamStation sleep therapy device the way you want with the connectivity that you want.



Ordering information

APAP/CPAP Sleep therapy

DreamStation Auto CPAP	Part number	With cellular modem	With WiFi modem
DreamStation Auto CPAP	DSX500S11	DSX500S11C05	DSX500S11W05
DreamStation Auto CPAP with humidifier	DSX500H11	DSX500H11C	DSX500H11W
DreamStation Auto CPAP with humidifier and heated tube	DSX500T11	DSX500T11C	DSX500T11W

DreamStation CPAP Pro	Part number	With cellular modem	With WiFi modem
DreamStation CPAP Pro	DSX400S11	DSX400S11C05	DSX400S11W05
DreamStation CPAP Pro with humidifier	DSX400H11	DSX400H11C	DSX400H11W
DreamStation CPAP Pro with humidifier and heated tube	DSX400T11	DSX400T11C	DSX400T11W

DreamStation CPAP	Part number	With cellular modem	With WiFi modem
DreamStation CPAP	DSX200S11	DSX200S11C05	DSX200S11W05
DreamStation CPAP with humidifier	DSX200H11	DSX200H11C05	DSX200H11W05
DreamStation CPAP with humidifier and heated tube	DSX200T11	DSX200T11C05	DSX200T11W05

Bi-level Sleep therapy

DreamStation Auto BiPAP	Part number	With cellular modem	With WiFi modem
DreamStation Auto BiPAP	DSX700S11	DSX700S11C05	DSX700S11W05
DreamStation Auto BiPAP with humidifier	DSX700H11	DSX700H11C	DSX700H11W
DreamStation Auto BiPAP with humidifier and heated tube	DSX700T11	DSX700T11C	DSX700T11W

DreamStation BiPAP Pro	Part number	With cellular modem	With WiFi modem
DreamStation BiPAP Pro	DSX600S11	DSX600S11C05	DSX600S11W05
DreamStation BiPAP Pro with humidifier	DSX600H11	DSX600H11C	DSX600H11W
DreamStation BiPAP Pro with humidifier and heated tube	DSX600T11	DSX600T11C	DSX600T11W

Ordering information (continued)

Sleep and Respiratory therapy

Humidification	Part number
DreamStation humidifier	DSXH
DreamStation water chamber	1122520
RP-DS Humidifier Dry Box Assembly	1120668
RP-DS Humidifier Flip Lid Seal	1120617
RP-DS Humidifier Dry Box Inlet Seal	1120613
Tubes	Part number
15mm Std Tube – DreamStation - RP	PR15
15mm Heated Tube	HT15
Filters	Part number
Pollen filter, reusable (1 per pack)	1122446
Ultra-fine filter, disposable (1 per pack)	1122447
Ultra-fine filter, disposable (2 per pack)	1122518
Ultra-fine filter, disposable (6 per pack)	1122519

Sleep and Respiratory therapy accessories

Power management	Part number
DreamStation 65W power supply*	1120136
DreamStation 80W power supply	1118499
Power cord	1005894
PAP Lithium Ion Battery Kit**	1124009
Li Ion battery, PAP device cable	1121736
Li Ion Battery, Cable Adapter (50 series adapter, 2/each)	1127799
Li Ion Battery Case	1126947

Sleep and Respiratory therapy accessories

Oximetry	Part number
Link module	1120293
Nonin SpO ₂ , assembly	1121694
Reuseable finger sensor (clip)	936
Oximetry finger sensor, adult, flex	953
Flexiwrap, sens tape, adult, 25/pk	954A
Reuseable finger sensor (soft), medium	1062050
Reuseable finger sensor (soft), small	1062051
Connectivity	Part number
Cellular modem, oximetry capable	100600C
Cellular modem, Non-oximetry	100605C
WiFi accessory, oximetry capable	100700W
WiFi accessory, Non-oximetry	100705W
Travel	Part number
Universal PAP/laptop travel briefcase	1114784
DreamStation replacement carrying case	1121162
Shielded DC cord	1120746
Shielded 12V DC power cord system (includes shielded DC cord and battery adapter cable)	1120747
DreamStation Travel case, 65W power supply* w/ USB port	1120135
DreamStation Travel case, 65W power supply* w/ USB port, 5/case	1124301
Mask and tubing bedside organizer	Part number
Bedside organizer	1122135
Bedside organizer, 10/case	1124224

*Used with flow generator only. Cannot be used with a humidifier.

**Includes a rechargeable lithium ion battery pack, a PAP device cable, and a battery case.

© 2016 Koninklijke Philips N.V. All rights reserved.
Specifications are subject to change without notice.



Caution: US federal law restricts these devices to sale by or on the order of a physician.

edoc KD 10/4/16 MCI 4107762
1010 Murry Ridge Lane, Murrysville, PA 15668
800-345-6443 • 724-387-4000

www.philips.com/respironics