



Philips InnoSpire Go  
Mask

**Adult large**

Soft seal mask  
For use with InnoSpire Go  
Replacement part

**HH1344**

## Aids in treatment

The detachable LiteTouch facemask has a soft contoured cushion that fits easily and rests gently on the face to provide a more comfortable treatment.\*

### **Soft seal technology**

- For adults who find it difficult using the mouthpiece

### **Rests gently on the face**

- With soft-seal technology\*\*

**PHILIPS**

Mask  
Adult large Soft seal mask, For use with InnoSpire Go, Replacement part

# Specifications

## Package contents

- Include: Large LiteTouch Mask

## Product details

- LiteTouch mask: For use with InnoSpire Go
- Packaging size: 82 mm x 97 mm x 100 mm (3.23" x

- 3.82" x 3.94")
- Packaged weight: 80 g
- Materials: Polyester, Silicone
- Size: Adult large
- Warranty: 1 year

HH1344/03

# Highlights

## Makes therapy easier



For adults who find it difficult using the mouthpiece, the LiteTouch masks make treatment easier.

## Rests gently on the face\*\*



The LiteTouch mask has a soft contoured cushion that rests gently on the face, so requires minimal pressure to apply.



Issue date 2018-07-26

Version: 6.0.2

12 NC: 0000 011 27875  
EAN: 87 12581 75432 7

© 2018 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\*\* Adapted from Erzinger et al. J Aerosol Med 2007 20 (Suppl1 S78-S84)

\*\* Quantifying facemask sealing efficiency when used on a valved holding chamber during simulated breathing. D von Hollen, E Lieberman, K Nikander, Philips Respironics, Respiratory Drug Delivery, Parsippany, NJ, USA. Presented at the Association of Asthma Educators (AAE) annual conference, 3 July