

# PHILIPS

Standard Mask


InnoSpire Go

Large LiteTouch Mask (5+ years)

Soft seal mask

For use with InnoSpire Go

HH1344



## Aids in treatment

The detachable LiteTouch facemask has a soft contoured cushion that fits easily and rests gently on the face to provide a more comfortable treatment.\*

### **Soft seal technology**

- For adults who find it difficult using the mouthpiece

### **Rests gently on the face**

- With soft-seal technology\*\*

# Highlights

## Makes therapy easier

For adults who find it difficult using the mouthpiece, the LiteTouch masks make treatment easier.

## Rests gently on the face\*\*

The LiteTouch mask has a soft contoured cushion that rests gently on the face, so requires minimal pressure to apply.

# Specifications

## Package contents

**Include:** Large LiteTouch Mask

## Product details

**Packaging size:** 82mm x 97mm x 100mm

**Packaged weight:** 80g

**LiteTouch mask:** For use with InnoSpire Go

**Materials:** Polyethylene terephthalate (PET), Silicone

**Size:** Adult large

**Warranty:** 1 year

## Material

**BPA:** Not manufactured with BPA



\*\* Adapted from Erzingerat al. J Aerosol Med 2007 20 (Suppl 1 S78-S84)

\*\*\* Quantifying facemask sealing efficiency when used on a valved holding chamber during simulated breathing. D von Hollen, E Lieberman, K Nikander, Philips Respironics, Respiratory Drug Delivery, Parsippany, NJ, USA. Presented at the Association of Asthma Educators (AAE) annual conference, July 3